

## OVERVIEW

### Protective Face Shields

Well suited for protection against COVID-19 and other pathogens. Ideal for healthcare and laboratory workers.

Made from PETG. Provides clarity of view. Visors are anti-fogging. Once size fits most. Includes a comfortable forehead cushion.

**CAPACITY: 10k Units Per Week.**

**MOQ: 10k**

**Price: \$4.89 per unit**

**Ships within 7-10 days of PO.**



## SPECIFICATIONS

**MATERIALS:** Shield: Food-grade PET with removable protective film  
Elastic: Polyester rubber  
Foam: High dense, open foam

**CHARACTERISTICS:**

- Anti-fog treatment
- Fiberglass-free
- Latex-free
- Lightweight Distortion-free
- Foam headband
- Optically clear

**DIMENSIONS:** Shield: 330 x 220 x 0.20 mm  
Foam: 260 x 30 x 30 mm  
Elastic: 340 x 150 mm

## CERTIFICATIONS

This product is manufactured in an FDA approved and CE certified factory.

Cleaning instructions: When disinfecting or sanitizing our face shields, the following products should be used: Steriplex SD, Safetec Surface Wipes, Perodox RTU, Prostat Sterile Wipes, hydrogen peroxide (pharmacy grade), bleach (household grade), or isopropyl alcohol (diluted with 30% water).

