

## OVERVIEW



### ASTM Level 1 – Surgical Mask

General Use Face Mask. COVID-19. Sterilized with ethylene oxide. Integrated layers work together to maximize breathability. Specially designed inner facing improves comfort. Soft comfortable ear loops provide easier donning. Offered in both blue & white, 20/Bag. ASTM Level 1.

**CAPACITY: 50k-100k Units Per Week.**

**MOQ: 50k**

**Price: \$1.65 per unit**

**Ships within 7-10 days of PO.**



## SPECIFICATIONS

**MATERIALS:**  
 Outer Layer: Leak proof non-woven fabric  
 Middle Layer: High density filler layer (melt spray nonwoven fabric)  
 Inner Layer: Spun-bonded non-woven fabric

**CHARACTERISTICS:**  
 Sterilized with ethylene oxide.  
 Soft comfortable ear loops provide easier donning.  
 Provides good bacterial and particle filtration.

≥95%	<b>BFE:</b> Bacterial filtration efficiency.
≥95%	<b>PFE:</b> Sub-micron particulates filtration efficient at 0.1 micron.
<4.0	<b>DP:</b> Differential pressure, mm H2O/cm2 (Breathability)
80 mm Hg	<b>RTP:</b> Resistance to penetration by synthetic blood, minimum pressure in mm Hg for pass results
Class 1	<b>FS:</b> Flame Spread

## CERTIFICATIONS

This product is manufactured in an FDA approved and CE certified factory.

\* These masks are intended for general, non-medical purposes, such as use in construction and other industrial applications. Because they are not intended for use in the diagnosis of disease or other conditions or in the cure, mitigation, treatment, or prevention of disease, FDA device marketing is not required and all the other requirements of the FD&C Act do not apply to manufacturers, importers and distributors of these products.

